



# Practical Pharmaceutical Analytical Chemistry - II

Second Level  
First Semester 2018-2019


Section 3

# Oxidation - Reduction Titration

## (Redox Titration)

# Redox Indicators

- ✓ They are highly colored organic compounds that change their color when being oxidized or reduced.
- ✓ The color change depends on the change in the redox potential (E) of the system during titration.
- ✓ The half-reaction responsible for the color change of the indicator can be written as follows:



➤ Requirements (specifications) of general redox indicators:

1. Has a standard redox potential **intermediate** between that of the sample and that of the titrant.

$$E^{\circ}_{\text{sample}} < E^{\circ}_{\text{indicator}} < E^{\circ}_{\text{titrant}} \text{ OR } E^{\circ}_{\text{sample}} > E^{\circ}_{\text{indicator}} > E^{\circ}_{\text{titrant}}$$

So that the titrant reacts first with the sample and then with the indicator at the end point.

2. Exhibits a **sharp, readily detectable** color change.

3. The **transition potential** of the indicator (i.e. the potential range at which the indicator changes its color) should be close to the potential at the equivalence point.

## ➤ Examples:

- 1,10-Phenanthroline indicator: { USP }

Used for the titration of ferrous {Fe<sup>2+</sup>} # cerric sulfate titrant {Ce(SO<sub>4</sub>)<sub>2</sub>}.

- Diphenylamine indicator:

Used for the titration of ferrous {Fe<sup>2+</sup>} # pot. dichromate titrant {K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub>}.

# Determination of Ferrous Salts ( $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$ )

**Experiment 1**


**Determination of Ferrous with  
0.1 N Ceric Sulphate using  
1,10-Phenanthroline Indicator**

# Principle


- ✓ Ferrous salts can be determined by titration with 0.1 N Ceric Sulphate using 1,10-Phenanthroline as a redox indicator.
- ✓ 1,10-Phenanthroline-ferrous complex (ferroin) is an intense red colored complex, which is reversibly oxidized (with strong oxidizing agents) to phenanthroline-ferric complex ion (ferrin), which has a pale blue color. The complex is used as an indicator in the titration of ferrous by ceric sulphate.



✓ During Titration:




✓ At End point:




# Procedure

**0.025 N St.  
ceric sulphate**




**Titrate Ferrous  
Sulphate with St.  
Ceric Sulphate**



**10 ml Ferrous sulphate sample  
+ 2 drops of 1,10 phenanthroline**

**At the E.P. the color in the conical  
flask changes from orange color to  
the first pale blue color**



# Calculations

## Calculation of equivalent:

$$\text{Each 1 ml 0.025 N Ce(SO}_4)_2 = \frac{278 \times 0.025}{1000} = 0.00695 \text{ g FeSO}_4 \cdot 7\text{H}_2\text{O.}$$

$$\text{Conc. of ferrous} = \frac{\text{E.p} \times \text{eq.} \times 100}{10} = \text{\% w/v}$$


## Experiment 2

**Determination of Ferrous with  
0.1 N Potassium Dichromate  
using Diphenylamine Indicator**

# Principle


✓ Ferrous salts can be determined by titration with 0.1 N Potassium Dichromate using Diphenylamine as a redox indicator.

✓ During Titration:



✓ At End point:



The diphenylamine (I) undergoes oxidation first into a colourless diphenylbenzidine (II) which is the real indicator and is reversibly further oxidized to diphenylbenzidine violet (III).



- ✓ The oxidation potential of the system diphenylamine/diphenylbenzidine (colourless)/diphenylbenzidine (violet) = 0.76 volt
- ✓ The oxidation potential of  $\text{Fe}^{3+}/\text{Fe}^{2+}$  = 0.77 volt
- ✓ The oxidation potential of  $\text{Cr}_2\text{O}_7^{2-}/2 \text{Cr}^{3+}$  = 1.36 volt
- It is obvious that there is an overlapping between  $E_{\text{ind}}$  &  $E_{\text{Fe}^{3+}/\text{Fe}^{2+}}$
- For this, diphenylamine is only able to function as indicator in this reaction when phosphoric acid is present in the solution.
- What is the role of Phosphoric acid?
- ❑ Formation of a colorless complex with the produced ferric ions leading to:
  1. Decreasing the molar concentration of Ferric and hence **reduces the actual potential of ( $\text{Fe}^{3+}/\text{Fe}^{2+}$ ) system** so that  $\text{Fe}^{2+}$  ion will be oxidized before the indicator.
  2. Removing the dark yellow color of  $\text{Fe}^{3+}$  ion giving a more clear color change.

# Electrochemical series ( $E^{\circ}$ values at $25^{\circ}\text{C}$ )


**TABLE 17.1 Standard Reduction Potentials at  $25^{\circ}\text{C}$**

	Reduction Half-Reaction	$E^{\circ}$ (V)			
<b>Stronger oxidizing agent</b> 	$\text{F}_2(\text{g}) + 2\text{e}^{-}$	$\longrightarrow 2\text{F}^{-}(\text{aq})$	2.87	<b>Weaker reducing agent</b> 	
	$\text{H}_2\text{O}_2(\text{aq}) + 2\text{H}^{+}(\text{aq}) + 2\text{e}^{-}$	$\longrightarrow 2\text{H}_2\text{O}(\text{l})$	1.78		
	$\text{MnO}_4^{-}(\text{aq}) + 8\text{H}^{+}(\text{aq}) + 5\text{e}^{-}$	$\longrightarrow \text{Mn}^{2+}(\text{aq}) + 4\text{H}_2\text{O}(\text{l})$	1.51		
	$\text{Cl}_2(\text{g}) + 2\text{e}^{-}$	$\longrightarrow 2\text{Cl}^{-}(\text{aq})$	1.36		
	$\text{Cr}_2\text{O}_7^{2-}(\text{aq}) + 14\text{H}^{+}(\text{aq}) + 6\text{e}^{-}$	$\longrightarrow 2\text{Cr}^{3+}(\text{aq}) + 7\text{H}_2\text{O}(\text{l})$	1.33		
	$\text{O}_2(\text{g}) + 4\text{H}^{+}(\text{aq}) + 4\text{e}^{-}$	$\longrightarrow 2\text{H}_2\text{O}(\text{l})$	1.23		
	$\text{Br}_2(\text{aq}) + 2\text{e}^{-}$	$\longrightarrow 2\text{Br}^{-}(\text{aq})$	1.09		
	$\text{Ag}^{+}(\text{aq}) + \text{e}^{-}$	$\longrightarrow \text{Ag}(\text{s})$	0.80		
	$\text{Fe}^{3+}(\text{aq}) + \text{e}^{-}$	$\longrightarrow \text{Fe}^{2+}(\text{aq})$	0.77		
	$\text{O}_2(\text{g}) + 2\text{H}^{+}(\text{aq}) + 2\text{e}^{-}$	$\longrightarrow \text{H}_2\text{O}_2(\text{aq})$	0.70		
	$\text{I}_2(\text{s}) + 2\text{e}^{-}$	$\longrightarrow 2\text{I}^{-}(\text{aq})$	0.54		
	$\text{O}_2(\text{g}) + 2\text{H}_2\text{O}(\text{l}) + 4\text{e}^{-}$	$\longrightarrow 4\text{OH}^{-}(\text{aq})$	0.40		
	$\text{Cu}^{2+}(\text{aq}) + 2\text{e}^{-}$	$\longrightarrow \text{Cu}(\text{s})$	0.34		
	$\text{Sn}^{4+}(\text{aq}) + 2\text{e}^{-}$	$\longrightarrow \text{Sn}^{2+}(\text{aq})$	0.15		
		<b><math>2\text{H}^{+}(\text{aq}) + 2\text{e}^{-}</math></b>	<b><math>\longrightarrow \text{H}_2(\text{g})</math></b>		<b>0</b>
		$\text{Pb}^{2+}(\text{aq}) + 2\text{e}^{-}$	$\longrightarrow \text{Pb}(\text{s})$		-0.13
	$\text{Ni}^{2+}(\text{aq}) + 2\text{e}^{-}$	$\longrightarrow \text{Ni}(\text{s})$	-0.26		
	$\text{Cd}^{2+}(\text{aq}) + 2\text{e}^{-}$	$\longrightarrow \text{Cd}(\text{s})$	-0.40		
	$\text{Fe}^{2+}(\text{aq}) + 2\text{e}^{-}$	$\longrightarrow \text{Fe}(\text{s})$	-0.45		
	$\text{Zn}^{2+}(\text{aq}) + 2\text{e}^{-}$	$\longrightarrow \text{Zn}(\text{s})$	-0.76		
	$2\text{H}_2\text{O}(\text{l}) + 2\text{e}^{-}$	$\longrightarrow \text{H}_2(\text{g}) + 2\text{OH}^{-}(\text{aq})$	-0.83		
	$\text{Al}^{3+}(\text{aq}) + 3\text{e}^{-}$	$\longrightarrow \text{Al}(\text{s})$	-1.66		
	$\text{Mg}^{2+}(\text{aq}) + 2\text{e}^{-}$	$\longrightarrow \text{Mg}(\text{s})$	-2.37		
	$\text{Na}^{+}(\text{aq}) + \text{e}^{-}$	$\longrightarrow \text{Na}(\text{s})$	-2.71		
<b>Weaker oxidizing agent</b>	$\text{Li}^{+}(\text{aq}) + \text{e}^{-}$	$\longrightarrow \text{Li}(\text{s})$	-3.04	<b>Stronger reducing agent</b>	




# Procedure

**0.1 N St.  
Potassium  
dichromate**



**Titrate Ferrous  
Sulphate with St.  
Potassium  
dichromate**



**10 ml Ferrous sulphate sample  
+ 5 ml dil. H<sub>2</sub>SO<sub>4</sub>  
+ 1 ml phosphoric acid  
+ 7 drops of diphenylamine**

The color in the conical flask changes from greenish blue color to greenish grey color then at the E.P. changes into intense violet blue color



# Calculations

## Calculation of equivalent:

$$\text{Each 1 ml 0.1 N K}_2\text{Cr}_2\text{O}_7 = \frac{278}{10 \times 1000} = 0.0278 \text{ g FeSO}_4 \cdot 7\text{H}_2\text{O}.$$

$$\text{Conc. of ferrous} = \frac{\text{E.P} \times \text{eq.} \times 100}{10} = \text{\% w/v}$$

*Thank You*

